
OPTIONS: for vertical autoclaves, need to be factory fitted.



Pt 100 Sensor. 2nd Probe for reading sample temperature.

To guarantee in all circumstances that the system has completed the cycle to the set time and temperature parameters, the Autester ST DRY PV II autoclave can be installed with a temperature probe that can be placed within the culture flask, residual packets etc. In this way proof that the programmed time and temperature for the run have been achieved.

When the probe is connected to the controller of the autoclave, it is automatically recognised by the microprocessor, whereby the real time and temperature values detected by the probe are displayed and printed. For autoclaves PRESOCLAVE II this Pt100 sensor only temperature reader.

Part No. 4001219

ACCESSORIES FOR VERTICAL AUTOCLAVES



Perforated AISI 316 stainless steel baskets.

Part No.	Ø cn	Height n	Suitable for autoclaves Part No.	Number of baskets
1004770	21.5	15	4001404	1
1000596	22.5	24	4047725	1
1000495	28	20	4002135 - 4002514	2
1000495	28	20	4001725 - 4002515	3
1000496	38	28	4002136 - 4002425	2
1000780	48	22	4002426	3



AISI 316 stainless steel wire dividers for autoclave baskets. Part No.

1001217 4 divider compartments for basket part no. 1000495.

1001218 6 divider compartments for basket part no. 1000496.

1001222 7 divider compartments for basket part no. 1000780.



Sterilisation Drums.

Part No.	Ø	Height cm	Suitable for autoclaves Part No.	Capacity Drums
1002418	28	18	4002135 - 4002514	2
1002418	28	18	4001725 - 4002515	3
1002419	35	27	4002136 - 4002425	2



Water distiller, specially for autoclaves "DEST-4".

Distillation capacity: 1.5 litres/hour. Internal capacity 4 litres. Case made from thermo resistant plastic.

Internal stainless steel. Cooling by forced air though a condenser. Without glass heater elements. Easy to use. Safety cut out if over temperature due to insufficient water. Comes complete with a 4 litre three way dispenser plastic bottle with filter. Note: for daily use it is recommended that the filter be changed every 3 months.

Part No. Conductivity Ø Height Weight Capacity Distillation litre/Hour 4001729 1.5 29 39 750 3.5

Spare: Dispenser spout with filter. Part No. 4001730



Detergent "BIO-SEL"

We recommend the use of BIO-Sel detergent for cleaning, it is also biodegradable and free from Phosphates.

Part No.

1000025 Capacity 1 litre. 1000026 Capacity 4 litres. 1000852 Box of 12 x 1 litre flasks.

1000853 Box of 4 x 4 litre flasks.

additional water required