


LABOREXIN OY

Lasinpuhaltamo Siljander

KYVETIT TEFLON TULPALLA - UV- VIS - NIR - IR

12,5x12,5x45mm


Kyveti tyyppi	Valotie mm	Tilavuus ml	Optinen lasi		UV-VIS		UV/VIS/NIR
			G	H	Q	I	
21	1	0.4	5121-G-1	5121-H-1	5121-Q-1	5121-I-1	
21	2	0.7	5121-G-2	5121-H-2	5121-Q-2	5121-I-2	
21	5	1.7	5121-G-5	5121-H-5	5121-Q-5	5121-I-5	
21	10	3.5	5121-G-10	5121-H-10	5121-Q-10	5121-I-10	
21	20	7.0	5121-G-20	5121-H-20	5121-Q-20	5121-I-20	
21	30	10.5	5121-G-30	5121-H-30	5121-Q-30	5121-I-30	
21	40	14.0	5121-G-40	5121-H-40	5121-Q-40	5121-I-40	
21	50	17.0	5121-G-50	5121-H-50	5121-Q-50	5121-I-50	



FLUOROMETRIKYVETIT TEFLON-TULPALLA

12,5x12,5x45mm

Kyveti tyyppi	Valotie mm	Tilavuus ml	Optinen lasi		UV-VIS		UV/VIS/NIR
			G	H	Q	I	
23	5	1.7	5123-G-5	5123-H-5	5123-Q-5	5123-I-5	
23	10	3.5	5123-G-10	5123-H-10	5123-Q-10	5123-I-10	
23	20	7.0	5123-G-20	5123-H-20	5123-Q-20	5123-I-20	
23	40	10.5	5123-G-40	5123-H-40	5123-Q-40	5123-I-40	
23	50	14.0	5123-G-50	5123-H-50	5123-Q-50	5123-I-50	




FLUOROMETRIKYVETTI ANAEROBINEN

28F	10	3.5	5128F-G-10	-	5128F-Q-10	-	
-----	----	-----	------------	---	------------	---	--

SEMI-MIKROKYVETIT TEFLON-TULPALLA

12,5x12,5x45mm

Kyveti tyyppi	Valotie mm	Tilavuus ml	Optinen lasi		UV-VIS		UV/VIS/NIR
			G	H	Q	I	
29	1	0.17	5129-G-1	5129-H-1	5129-Q-1	5129-I-1	
29	2	0.35	5129-G-2	5129-H-2	5129-Q-2	5129-I-2	
29	5	0.7	5129-G-5	5129-H-5	5129-Q-5	5129-I-5	
29	10	1.4	5129-G-10	5129-H-10	5129-Q-10	5129-I-10	
29	20	2.8	5129-G-20	5129-H-20	5129-Q-20	5129-I-20	
29	30	4.2	5129-G-30	5129-H-30	5129-Q-30	5129-I-30	
29	40	5.6	5129-G-40	5129-H-40	5129-Q-40	5129-I-40	
29	50	7.0	5129-G-50	5129-H-50	5129-Q-50	5129-I-50	



SEMI-MIKROKYVETIT TEFLON-TULPALLA, "SELF MASKING"

29B	5	0.7	5129B-G-5	5129B-H-5	5129B-Q-5	5129B-I-5	
29B	10	1.4	5129B-G-10	5129B-H-10	5129B-Q-10	5129B-I-10	
29B	20	2.8	5129B-G-20	5129B-H-20	5129B-Q-20	5129B-I-20	
29B	40	5.6	5129B-G-40	5129B-H-40	5129B-Q-40	5129B-I-40	



SEMI-MIKRO-FLUOROMETRIKYVETTI TEFLON-TULPALLA

29F	10	1.4	5129F-G-10	-	5129F-Q-10	-	
-----	----	-----	------------	---	------------	---	--

Puh. 09-780 633, Fax. 09-781 393

E-mail: laborexin@laborexin.fi www.laborexin.fi