


# LABOREXIN OY

## Lasinpuhaltamo Siljander

### MIKRO-KYVETIT TEFLON TULPALLA

12,5x12,5x45mm, UV - VIS - NIR - IR


Kyveti tyyppi	Valotie mm	Tilavuus ml	Optinen lasi		UV-VIS		UV/VIS/NIR
			G	H	Q	I	
30	1	0.8	5130-G-1	5130-H-1	5130-Q-1	5130-I-1	
30	2	0.17	5130-G-2	5130-H-2	5130-Q-2	5130-I-2	
30	5	0.35	5130-G-5	5130-H-5	5130-Q-5	5130-I-5	
30	10	0.7	5130-G-10	5130-H-10	5130-Q-10	5130-I-10	
30	20	1.4	5130-G-20	5130-H-20	5130-Q-20	5130-I-20	
30	30	2.1	5130-G-30	5130-H-30	5130-Q-30	5130-I-30	
30	40	2.8	5130-G-40	5130-H-40	5130-Q-40	5130-I-40	
30	50	3.5	5130-G-50	5130-H-50	5130-Q-50	5130-I-50	



### MIKRO-KYVETIT TEFLON TULPALLA, "SELF MASKIN"

12,5x12,5x45mm, UV - VIS - NIR - IR

Kyveti tyyppi	Valotie mm	Tilavuus ml	Optinen lasi		UV-VIS		UV/VIS/NIR
			G	H	Q	I	
30B	5	0.35	5130B-G-5	5130B-H-5	5130B-Q-5	5130B-I-5	
30B	10	0.7	5130B-G-10	5130B-H-10	5130B-Q-10	5130B-I-10	
30B	20	1.4	5130B-G-20	5130B-H-20	5130B-Q-20	5130B-I-20	
30B	40	2.8	5130B-G-40	5130B-H-40	5130B-Q-40	5130B-I-40	




### MIKRO-FLUOROMETRIKYVETTI TEFLON TULPALLA

30F	10	0.7	5130F-G-10	-	5130F-Q-10	-	
-----	----	-----	------------	---	------------	---	--

### LÄPIVIRTAUSKYVETIT MIKRO- UV - VIS - NIR - IR

12,5x12,5x45mm

Kyveti tyyppi	Valotie mm	Tilavuus ml	Optinen lasi		UV-VIS		UV/VIS/NIR
			G	H	Q	I	
44	5	0.9	5144-G-5	5144-H-5	5144-Q-5	5144-I-5	
44	10	1.8	5144-G-10	5144-H-10	5144-Q-10	5144-I-10	
44	20	3.8	5144-G-20	5144-H-20	5144-Q-20	5144-I-20	
44	40	7.2	5144-G-40	5144-H-40	5144-Q-40	5144-I-40	




### LÄPIVIRTAUSKYVETTI MIKRO- UV - VIS - NIR - IR, "SELF MASKING"

44B	10	1.8	5144B-G-10	5144-H-10	5144-Q-10	5144-I-10	
-----	----	-----	------------	-----------	-----------	-----------	--

### KYVETIT VAIPALLA, UV - VIS - NIR - IR

12,5x12,5x45mm

Kyveti tyyppi	Valotie mm	Tilavuus ml	Optinen lasi		UV-VIS		UV/VIS/NIR
			G	H	Q	I	
61	10	2.6	5161-G-10	5161-H-10	5161-Q-10	5161-I-10	
61	20	5.2	5161-G-20	5161-H-20	5161-Q-20	5161-I-20	
61	40	10.4	5161-G-40	5161-H-40	5161-Q-40	5161-I-40	



Holkkitie 14, PL 224, 00811 Helsinki

Puh. 09-780 633, Fax. 09-781 393

E-mail: laborexin@laborexin.fi www.laborexin.fi