## COMECTA: Microscope biological inverted, trinocular "XDS-1B"



## **APPLICATIONS**

Optical microscope with many benefits, specially designed for biomedicine, biology and for culture analysis in matrix flasks and in microtitre plates, flasks, culture vessels, tubes etc. through the use of its objectives and long working distance condenser, technique uses brightfield or phase contrast.

## **FEATURES**

Robust stable design with lateral arm wrist supports for maximum comfort and easy manipulation.

Trinocular head, inclined at 45°.

Interpupillary distance adjustable from 55 to 77 mm.

**Diopter correction** of  $\pm 5$ dp in the left eyepiece, used to correct differences in sharpness.

**Pair of DIN** wide field **WF 10x** (20 Ø), pair of WF 16x (14 Ø), plan eyepiece PF 10x (17.4 Ø). **Reticule** with 0.1 mm/div, telescopic central eyepiece phase contrast DZ11x DZ11x and trinocular micrograph with photography eyepiece **S 5x** (9.5 Ø).

**Objectives DIN planachromatic brightfield,long working distance LWD PL** of 10x (A.N. 0.25/working distance W.D.=7.9 mm), **25x** (A.N. 0.40/W.D.=5 mm), **40x** (A.N. 0.40/W.D.=3 mm) and on **negative phase contrast LWD PL ph of 25x** (A.N. 0.25/W.D.=5 mm), colour coded for rapid identification. Maximum thickness acceptable for receptacle bases 1.5 mm.

Quadruple objectives turret with rotation in both directions with a roller baring and notch stop. Rack and pinion height adjustment mechanism with bilateral coaxial focus controls: micrometer of 0.2 mm rotation graduated at 0.002 mm increments and macrometer with adjustable tension.

Total working distance 12 mm.

**Condenser with long working distance (55.2 mm) LWD 0.3 A.N.**, brightfield and phase contrast. Precentred, with iris diaphragm and flip out phase contrast. **Including a phase ring of 25x.** Rack and pinions height mechanism or lateral removal for receptacles of heights up to 145 mm.

Large stage: dimensions 180x155 mm. Includes 2 holding clips for tubes and slides and 4 insert discs of varying size for optimum accommodation of differing sample vessels: plates and or culture flasks, microtitre plates, etc.

Mechanical sub-stage with graduated scale of 1 mm and vernier of 0.1 mm, made up of a smooth mechanism (x-y) rack and pinion drive of 50 mm longitudinal and 75 mm transversal.

**Illumination system, diascopic** transmitted light by a halogen lamp of 6 V / 20W, precentred by external controls and controllable intensity.

The unit includes a field diaphragm, filter carousel of green, yellow, blue and grey.

MODEL	Part No.	Height / Width / Depth exterior cm	Weight Kg
XDS-1B	5901998	50 28 30	15



Comes complete with 2 flasks, 2 culture plates, wooden holder for eyepieces, objectives and accessories, 2 lamps, 2 spare fuses and a protective cover.

ACCESSORIES Objectives planachromatic lo Negative phase contrast: Objective LWD PL ph 10x (A.N. 0.25 160/1.5 W.D.=7.9 Objective LWD PL ph 40) (A.N. 0.60 160/1.5 W.D.=3 n	mm) Part No. <mark>5313190</mark> (	5313190	5313192
Phase contrast ring 10x.	Fait NU. 33 13 192	5313191	5313193
	Part No. <b>5313191</b>		
Phase contrast ring 40x.	Part No. <b>5313193</b>		
<b>CCD video camera eyepie</b> S0.65x.			A
<b>Adapter thread C</b> for video camera CCD.	Part No. <b>5313220</b> Part No. <b>5313219</b>		men
Video camera CCD.	Part No. <b>5435038</b>	5313219	5435038
Motorised camera back Pentax de 35 mm, mount T2 and electronic trigger cable.		PEN PEN	TAX
Photographic eyepiece Si	Part No. 5313216		
<b>SPARES</b> Halogen lamp 6 V / 20 W.	Part No. <b>5313023</b>	5313216	Dare-

258